


Anti-C-Abl Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P00519

Item Number: CM0017688



Introduction

Rabbit polyclonal antibody against human c-Abl (ABL1), validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat. Ideal for studying cell growth, survival, DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	c-Abl Rabbit pAb
Remarks/Alias	ABL; JTK7; p150; c-ABL; v-abl; CHD5KM; c-ABL1; BCR-ABL; bcr/abl; c-Abl
Species	Human
GeneID (Human)	25
GeneID	25
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1000-1130 of human c-Abl (NP_005148.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00519
Protein Weight	123kDa
Shipping	Ice bag