


Anti-C-Fos Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P01100

Item Number: CM0007000



Introduction

Recombinant rabbit monoclonal c-Fos antibody validated for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Reacts with human, mouse, and rat c-Fos. Suitable for studying AP-1 transcription factor signaling.

[Learn More](#)

Parameter	Value
Product Name	c-Fos Rabbit mAb
Remarks/Alias	p55; AP-1; C-FOS; c-Fos
Species	Human
Gene ID (Human)	2353
Gene ID	2353
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 211-281 of human c-Fos (NP_005243.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01100
Protein Weight	41kDa
Shipping	Ice bag