


Anti-C-Jun Rabbit Polyclonal Antibody For Wb And Elisa - P05412

Item Number: CM0018168



Introduction

c-Jun Rabbit Polyclonal Antibody (pAb) is validated for Western blot (WB) and ELISA applications. It detects c-Jun transcription factor in human, mouse, and rat samples, making it suitable for AP-1 pathway and cancer research.

[Learn More](#)

Parameter	Value
Product Name	c-Jun Rabbit pAb
Remarks/Alias	AP1; p39; AP-1; cJUN; c-Jun
Species	Human
GeneID (Human)	3725
GeneID	3725
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-70 of human c-Jun (NP_002219.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05412
Protein Weight	36kDa
Shipping	Ice bag