


Anti-C-Jun Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P05412

Item Number: CM0018659



Introduction

c-Jun Rabbit Polyclonal Antibody (pAb) validated for WB, IF/ICC, ELISA. Detects human and mouse c-Jun (36kDa), a transcription factor involved in AP-1 signaling and oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	c-Jun Rabbit pAb
Remarks/Alias	AP1; p39; AP-1; cJUN; c-Jun
Species	Human
Gene ID (Human)	3725
Gene ID	3725
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-230 of c-Jun (NP_002219.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05412
Protein Weight	36kDa
Shipping	Ice bag