


Anti-C-Myc Polyclonal Antibody For Wb, Elisa - P01106

Item Number: CM0022937



Introduction

Rabbit polyclonal antibody against human c-Myc (MYC, bHLHe39), suitable for Western blot and ELISA applications. Detects human and mouse c-Myc protein (~51 kDa). Useful in cancer, stem cell, and transcription research.

[Learn More](#)

Parameter	Value
Product Name	c-Myc Rabbit pAb
Remarks/Alias	MRTL; MYCC; c-Myc; bHLHe39
Species	Human
GeneID	4609
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 347-439 of human c-Myc (NP_001341799.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01106
Protein Weight	51kDa
Shipping	Ice bag