


Anti-C-Myc Rabbit Polyclonal Antibody For Wb, Elisa - P01106

Item Number: CM0016390



Introduction

Rabbit polyclonal antibody against c-Myc transcription factor, verified for WB and ELISA in human and mouse. Suitable for oncogene signaling, transcription regulation, and stem cell research.

[Learn More](#)

Parameter	Value
Product Name	c-Myc Rabbit pAb
Remarks/Alias	MRTL; MYCC; c-Myc; bHLHe39
Species	Human
GeneID (Human)	4609
GeneID	4609
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 347-439 of human c-Myc (NP_001341799.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01106
Protein Weight	51kDa
Shipping	Ice bag