

Anti-C-Peptide Rabbit Monoclonal Antibody For If/Icc, If-P, Ihc-P, Elisa - P01308

Item Number: CM0005156



Introduction

Rabbit monoclonal antibody targeting C-Peptide, validated for IF/ICC, IF-P, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for diabetes research and monitoring beta-cell function.

[Learn More](#)

Parameter	Value
Product Name	C-Peptide Rabbit mAb
Remarks/Alias	IDDM; ILPR; IRDN; IDDM1; IDDM2; PNDM4; MODY10; C-Peptide
Species	Human
Gene ID (Human)	3630
Gene ID	3630
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human INS (NP_000198.1)
Source	Rabbit
Category	Monoclonal Antibody
Application	IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01308
Protein Weight	12kDa
Shipping	Ice bag