


Anti-C-Reactive Protein (Crp) Monoclonal Antibody For Wb, Ihc-P, Elisa - P02741

Item Number: CM0003850



Introduction

Rabbit monoclonal antibody against C-reactive protein (CRP), raised with a synthetic peptide. Validated for WB, IHC-P, and ELISA applications, cross-reacts with human and mouse CRP. Ideal for inflammation and immune response research.

[Learn More](#)

Product Name	C-Reactive Protein (CRP) Rabbit mAb
Remarks/Alias	PTX1; C-Reactive Protein (CRP)
Species	Human
Gene ID (Human)	1401
Gene ID	1401
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 125-224 of human C-Reactive Protein (CRP) (P02741).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02741
Protein Weight	23kDa
Shipping	Ice bag