


Anti-C-Reactive Protein (Crp) Polyclonal Antibody For Wb, Elisa - P02741

Item Number: CM0014357



Introduction

Rabbit anti-human C-reactive protein (CRP) polyclonal antibody validated for Western blot (WB) and ELISA, with mouse cross-reactivity. Ideal for acute-phase protein and host defense research and IVD assay development.

[Learn More](#)

Parameter	Specification
Product Name	C-Reactive Protein (CRP) Rabbit pAb
Remarks/Alias	PTX1; C-Reactive Protein (CRP)
Species	Human
Gene ID (Human)	1401
Gene ID	1401
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-50 of human C-Reactive Protein (CRP) (NP_000558.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P02741
Protein Weight	23kDa
Shipping	Ice bag