


Anti-C-Reactive Protein (Crp) Polyclonal Antibody For Wb, If/Icc, Elisa - P02741

Item Number: CM0014943



Introduction

C-Reactive Protein (CRP) rabbit polyclonal antibody for Western blot (WB), immunofluorescence/immunocytochemistry (IF/ICC), and ELISA. Cross-reacts with human, mouse, and rat. Ideal for acute phase and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	C-Reactive Protein (CRP) Rabbit pAb
Remarks/Alias	PTX1; C-Reactive Protein (CRP)
Species	Human
GeneID (Human)	1401
GeneID	1401
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human C-Reactive Protein (CRP) (NP_000558.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02741
Protein Weight	23kDa
Shipping	Ice bag