


Anti-C10Orf2 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Rr1

Item Number: CM0028968



Introduction

High-quality anti-C10orf2 rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, rat Twinkle mtDNA helicase, useful for mitochondrial DNA replication and repair studies.

[Learn More](#)

Parameter	Value
Product Name	C10orf2 Rabbit pAb
Remarks/Alias	PEO; PEO1; SCA8; ATXN8; IOSCA; PEOA3; SANDO; TWINL; MTDPS7; PRLT55; C10orf2
Species	Human
Gene ID (Human)	56652
Gene ID	56652
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 135-335 of human C10orf2 (NP_068602.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96RR1
Protein Weight	77kDa
Shipping	Ice bag