

Anti-C10Orf32 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96B45

Item Number: CM0029322

Introduction

Rabbit polyclonal antibody against human C10orf32 (BORCS7), validated for WB, IF/ICC, and ELISA. Cross-reacts with rat. Targets a lysosome-associated protein involved in intracellular trafficking and lysosome positioning.

[Learn More](#)

Parameter	Value
Product Name	C10orf32 Rabbit pAb
Remarks/Alias	C10orf32
Species	Human
GeneID	119032
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human C10orf32 (NP_653192.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96B45
Protein Weight	12kDa
Shipping	Ice bag