

Anti-C10Orf67 Rabbit Polyclonal Antibody For Wb, Elisa - Q8lyj2

Item Number: CM0021990

Introduction

Anti-C10orf67 rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Targets human C10orf67, cross-reacts with mouse and rat. Suitable for mitochondrial protein research.

[Learn More](#)

Parameter	Value
Product Name	C10orf67 Rabbit pAb
Remarks/Alias	C10orf115; LINC01552; C10orf67
Species	Human
GeneID (Human)	256815
GeneID	256815
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 68-243 of C10orf67 (NP_001338235.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8LYJ2
Protein Weight	64kDa
Shipping	Ice bag