

Anti-C12Orf10 Rabbit Polyclonal Antibody For IHC-P And Elisa - Q9Hb07

Item Number: CM0018865

Introduction

C12orf10 Rabbit Polyclonal Antibody for IHC-P and ELISA, reacting with human, mouse, and rat MYG1 exonuclease. Suitable for investigating RNA processing and nucleo-mitochondrial crosstalk.

[Learn More](#)

Parameter	Value
Product Name	C12orf10 Rabbit pAb
Remarks/Alias	MYG; Gamm1; MST024; MSTP024; C12orf10
Species	Human
Gene ID (Human)	60314
Gene ID	60314
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 48-376 of human C12orf10 (NP_067653.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HB07
Protein Weight	42kDa
Shipping	Ice bag