

Anti-C14Orf166/Rtraf Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y224

Item Number: CM0027758



Introduction

High-quality rabbit polyclonal antibody targeting human C14orf166 (RTRAF), validated for WB and ELISA. Reacts with human and mouse samples. Protein weight 28 kDa. Ideal for studying tRNA splicing and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	C14orf166 Rabbit pAb
Remarks / Alias	CLE; CLE7; hCLE; CGI99; RLLM1; hCLE1; CGI-99; LCRP369; C14orf166
Species	Human
GeneID	51637
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-244 of human C14orf166 (NP_057123.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y224
Protein Weight	28kDa
Shipping	Ice bag