


Anti-C17Orf70 Rabbit Polyclonal Antibody For Wb And Elisa - Q0Vg06

Item Number: CM0030575



Introduction

Anti-C17orf70 rabbit polyclonal antibody validated for WB and ELISA, recognizing human, mouse, and rat C17orf70 (FAAP100). Useful for studying DNA damage repair and Fanconi anemia pathways.

[Learn More](#)

Parameter	Value
Product Name	C17orf70 Rabbit pAb
Remarks/Alias	C17orf70
Species	Human
GeneID	80233
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 652-881 of human C17orf70 (NP_079437.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q0VG06
Protein Weight	93kDa
Shipping	Ice bag