

Anti-C18Orf8 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Dm3

Item Number: CM0020235



Introduction

Rabbit polyclonal antibody targeting C18orf8/RMC1 (MIC1), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. For autophagy and endosomal trafficking studies.

[Learn More](#)

Parameter	Details
Product Name	C18orf8 Rabbit pAb
Remarks/Alias	MIC1; Mic-1; WDR98; C18orf8; HsT2591
Species	Human
Gene ID (Human)	29919
Gene ID	29919
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human C18orf8 (NP_037458.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96DM3
Protein Weight	75kDa
Shipping	Ice bag