


Anti-C19Orf10 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q969H8

Item Number: CM0017073



Introduction

C19orf10 Rabbit Polyclonal Antibody (MYDGF) is designed for WB, IHC-P, ELISA. Recognizes human, mouse, rat. Target is myeloid-derived growth factor, involved in cardiac protection and angiogenesis. SWISS: Q969H8.

[Learn More](#)

Parameter	Value
Product Name	C19orf10 Rabbit pAb
Remarks/Alias	IL25, IL27, SF20, IL27w, C19orf10, R33729_1, EUROIMAGE1875335
Species	Human
GeneID	56005
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-173 of human C19orf10 (NP_061980.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969H8
Protein Weight	19kDa
Shipping	Ice bag