

Anti-C1Galt1C1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Eu7

Item Number: CM0028002

Introduction

C1GALT1C1 Rabbit pAb is a polyclonal antibody validated for WB and ELISA, targeting human C1GALT1C1 (COSMC) chaperone involved in O-glycan biosynthesis. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	C1GALT1C1 Rabbit pAb
Remarks/Alias	TNPS; COSMC; MST143; C1GALT2; HSPC067; C38H2-L1; C1GALT1C1
Species	Human
Gene ID (Human)	29071
Gene ID	29071
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96EU7
Protein Weight	36kDa
Shipping	Ice bag