

Anti-C1Orf146/Spo16 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q3Kqp7

Item Number: CM0023768

Introduction

Rabbit polyclonal antibody targeting mouse SPO16/C1orf146, validated for WB, IF-P, and ELISA. Detects mouse and rat SPO16 (~21 kDa). Ideal for meiosis and synaptonemal complex research.

[Learn More](#)

Parameter	Value
Product Name	C1orf146 Rabbit pAb
Remarks/Alias	SCRE, Spo16, C1orf146
Species	Mouse
Gene ID	76421
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 47-173 of mouse C1orf146 (NP_877422.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q3KQP7
Protein Weight	21kDa
Shipping	Ice bag