


Anti-C1Qa Polyclonal Antibody For Wb, Elisa - P02745

Item Number: CM0021099



Introduction

Rabbit polyclonal antibody against human C1QA for WB and ELISA. Validated in human, mouse, rat. Recognizes complement C1q subcomponent subunit A, a key initiator of the classical complement pathway.

[Learn More](#)

Parameter	Value
Product Name	C1QA Rabbit pAb
Remarks/Alias	C1QD1; C1QA
Species	Human
Gene ID (Human)	712
Gene ID	712
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-245 of human C1QA (NP_001334394.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02745
Protein Weight	26kDa
Shipping	Ice bag