

Anti-C1Qa Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P98086

Item Number: CM0023445

Introduction

Rabbit polyclonal antibody targeting mouse C1qa, validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat C1qa (26kDa). Ideal for complement classical pathway and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	C1qa Rabbit pAb
Remarks/Alias	C1q; Adic; C1qa
Species	Mouse
Gene ID	12259
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-245 of mouse C1qa (NP_031598.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98086
Protein Weight	26kDa
Shipping	Ice bag