


Anti-C1Qb Polyclonal Antibody For Wb, Elisa - P02746

Item Number: CM0030836



Introduction

Rabbit polyclonal antibody against human C1QB, validated for WB and ELISA with cross-reactivity to human and mouse. Ideal for complement classical pathway research and immune complex studies.

[Learn More](#)

Parameter	Value
Product Name	C1QB Rabbit pAb
Remarks/Alias	C1QD2; C1QB
Species	Human
GeneID (Human)	713
GeneID	713
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 104-253 of human C1QB (NP_000482.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02746
Protein Weight	27kDa
Shipping	Ice bag