


Anti-C1Qc Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P02747

Item Number: CM0026324



Introduction

Rabbit polyclonal antibody targeting human C1QC (complement C1q subcomponent subunit C). Suitable for WB, IHC-P, ELISA. Essential for classical complement pathway and innate immunity studies. Immunogen: synthetic peptide aa 50-150. Cross-reacts with human. MW: 26kDa.

[Learn More](#)

Parameter	Value
Product Name	C1QC Rabbit pAb
Remarks/Alias	C1QG; C1Q-C; C1QD3; C1QC
Species	Human
GeneID (Human)	714
GeneID	714
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human C1QC (NP_001107573.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P02747
Protein Weight	26kDa
Shipping	Ice bag