

Anti-C1Qtnf5 Polyclonal Antibody For Wb And Elisa - Q9Bxj0

Item Number: CM0029439



Introduction

Rabbit polyclonal antibody against human C1QTNF5 (CTRP5). Suitable for WB and ELISA, cross-reacts with mouse and rat. Recombinant protein immunogen. For research in extracellular matrix and retinal degeneration.

[Learn More](#)

Parameter	Value
Product Name	C1QTNF5 Rabbit pAb
Remarks/Alias	MFRP, CTRP5, C1QTNF5
Species	Human
Gene ID (Human)	114902
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-243 of human C1QTNF5 (NP_056460.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXJ0
Protein Weight	25kDa
Shipping	Ice bag