

Anti-C1R Polyclonal Antibody For If/Icc And Elisa - P00736

Item Number: CM0025934



Introduction

Rabbit polyclonal antibody against human C1R (Complement C1r subcomponent, EDS8, EDSPD1). Suitable for IF/ICC and ELISA, cross-reacts with human and mouse. Research use for complement pathway and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	C1R Rabbit pAb
Remarks/Alias	EDS8; EDSPD1; C1R
Species	Human
GeneID (Human)	715
GeneID	715
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 466-705 of human C1R (NP_001724.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00736
Protein Weight	80kDa
Shipping	Ice bag