

Anti-C1R Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - P00736

Item Number: CM0019982



Introduction

Rabbit polyclonal antibody against human C1R, validated for IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for complement pathway research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	C1R Rabbit pAb
Remarks/Alias	EDS8; EDSPD1
Species	Human
GeneID (human)	715
GeneID	715
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 18-190 of human C1R(NP_001724.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00736
Protein Weight	80kDa
Shipping	Ice bag