

Anti-C1S Polyclonal Antibody For Wb, Elisa - P09871

Item Number: CM0028101



Introduction

Rabbit polyclonal antibody against human C1S, suitable for WB and ELISA, with cross-reactivity to rat. Useful for classical complement pathway research. Protein weight: 77kDa.

[Learn More](#)

Parameter	Value
Product Name	C1S Rabbit pAb
Remarks/Alias	EDSPD2; C1S
Species	Human
GeneID (Human)	716
GeneID	716
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 479-688 of human C1S (NP_958850.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P09871
Protein Weight	77kDa
Shipping	Ice bag