

Anti-C21Orf33 Polyclonal Antibody For Wb, If/Icc, Elisa - P30042

Item Number: CM0026809

Introduction

High-quality rabbit polyclonal antibody targeting human C21orf33, validated for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat samples. Suitable for protein detection and IVD raw material sourcing.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] C21orf33 Rabbit pAb
Remarks/Alias	ES1; HES1; KNPH; KNPI; GT335; GATD3A; C21orf33; 33
Species	Human
GeneID (Human)	8209
GeneID	8209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-268 of human C21orf33 (NP_004640.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30042
Protein Weight	24kDa/28kDa
Shipping	Ice bag