

Anti-C21Orf33 Rabbit Polyclonal Antibody (Ko Validated) For Wb, Elisa - P30042

Item Number: CM0021861

Introduction

KO validated rabbit polyclonal antibody targeting human C21orf33 (ES1, HES1, KNPH). Suitable for WB and ELISA; cross-reacts with mouse and rat. Ideal for exploring C21orf33 function in cellular contexts.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] C21orf33 Rabbit pAb
Remarks/Alias	ES1; HES1; KNPH; KNPI; GT335; GATD3A; C21orf33; 33
Species	Human
Gene ID (Human)	8209
Gene ID	8209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-268 of human C21orf33 (NP_004640.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30042
Protein Weight	24kDa/28kDa
Shipping	Ice bag