

Anti-C2Cd5 Rabbit Polyclonal Antibody For Wb And Elisa - Q86Ys7

Item Number: CM0016149

Introduction

Rabbit polyclonal antibody targeting human C2CD5, suitable for western blot and ELISA. Cross-reacts with mouse and rat. Ideal for glucose transport and insulin signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	C2CD5 Rabbit pAb
Remarks/Alias	CDP138; KIAA0528; C2CD5
Species	Human
Gene ID (Human)	9847
Gene ID	9847
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 732-1051 of human C2CD5 (NP_001273103.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86YS7
Protein Weight	110kDa
Shipping	Ice bag