

Anti-C3 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P01024

Item Number: CM0012812



Introduction

Rabbit monoclonal antibody against human C3, targeting the C3c region (aa 1364-1663). Suitable for WB, IF/ICC, ELISA and cross-reacts with mouse and rat. Ideal for complement pathway and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	C3 Rabbit mAb
Remarks/Alias	ASP, C3a, C3b, AHUS5, ARMD9, CPAMD1, HEL-S-62p
Species	Human
GeneID (Human)	718
GeneID	718
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1364-1663 of human C3c (NP_000055.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01024
Protein Weight	187kDa
Shipping	Ice bag