


Anti-C3 Rabbit Polyclonal Antibody For Wb, Elisa - P01027

Item Number: CM0021758



Introduction

C3 Rabbit Polyclonal Antibody (pAb) validated for WB and ELISA. Detects mouse and rat Complement C3 (186 kDa) with high specificity. Ideal for complement pathway research and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	C3 Rabbit pAb
Remarks/Alias	ASP; Plp; HSE-MSF; C3
Species	Mouse
Gene ID (Human)	718
Gene ID	12266
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1159-1286 of mouse C3 (NP_033908.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P01027
Protein Weight	186kDa
Shipping	Ice bag