

Anti-C3C Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P01024

Item Number: CM0023315



Introduction

C3c Rabbit pAb is a polyclonal antibody raised against human C3c for WB, IF-P, and ELISA applications. It targets the complement C3 protein, a key component of innate immunity.

[Learn More](#)

Parameter	Value
Product Name	C3c Rabbit pAb
Remarks/Alias	ASP; C3a; C3b; AHUS5; ARMD9; CPAMD1; HEL-S-62p; C3c
Species	Human
GeneID (Human)	718
GeneID	718
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1364-1663 of human C3c. (NP_000055.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	P01024
Protein Weight	187kDa
Shipping	Ice bag