

Anti-C3Orf38 Rabbit Monoclonal Antibody For Wb And Elisa - Q5Jpi3

Item Number: CM0002123

Introduction

High-quality C3orf38 rabbit monoclonal antibody for human apoptosis research. Validated for Western blot and ELISA, detects uncharacterized protein C3orf38 (38kDa, SWISS: Q5JPI3).

[Learn More](#)

Parameter	Value
Product Name	C3orf38 Rabbit mAb
Species	Human
Gene ID (Human)	285237
Gene ID	285237
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 51-200 of human C3orf38 (Q5JPI3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5JPI3
Protein Weight	38kDa
Shipping	Ice bag