

Anti-C4Bpa Polyclonal Antibody For Wb And Elisa - P04003

Item Number: CM0026339



Introduction

High-quality rabbit polyclonal antibody targeting C4b-binding protein alpha chain (C4BPA), suitable for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Ideal for complement pathway research.

[Learn More](#)

Parameter	Value
Product Name	C4BPA Rabbit pAb
Remarks/Alias	PRP; C4BP; C4BPA
Species	Human
GeneID (human)	722
GeneID	722
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 500 to the C-terminus of human C4BPA (NP_000706.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04003
Protein Weight	67kDa
Shipping	Ice bag