

Anti-C4Bpb Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P20851

Item Number: CM0025811

Introduction

Rabbit polyclonal antibody against human C4BPB, suitable for WB, IF-P, and ELISA. Cross-reacts with mouse and rat. Targets the C4b-binding protein beta chain, a complement pathway regulator.

[Learn More](#)

Parameter	Value
Product Name	C4BPB Rabbit pAb
Remarks/Alias	C4BP; C4BPB
Species	Human
GeneID	725
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-251 of human C4BPB (NP_001017364.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20851
Protein Weight	28kDa
Shipping	Ice bag