

Anti-C4C Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P0C0L5

Item Number: CM0024145



Introduction

Rabbit polyclonal antibody targeting human C4c (Complement C4-B), validated for WB, IF/ICC, IF-P, ELISA. Cross-reacts with mouse and rat. Essential for complement pathway and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	C4c Rabbit pAb
Remarks/Alias	CH; C4F; CO4; C4B1; C4B2; C4B3; C4B5; C4BD; C4B12; C4B_2; CPAMD3; C4c
Species	Human
GeneID (Human)	721
GeneID	721
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1337-1446 of human C4c (NP_001002029.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0C0L5
Protein Weight	192kDa
Shipping	Ice bag