

Anti-C5Ar2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9P296

Item Number: CM0016102

Introduction

C5AR2 Rabbit Polyclonal Antibody (C5L2/GPR77) for WB, ELISA applications. Recognizes human and mouse C5a receptor 2, a GPCR mediating anaphylatoxin signaling. Suitable for inflammation and immunology studies.

[Learn More](#)

Parameter	Value
Product Name	C5AR2 Rabbit pAb
Remarks/Alias	C5L2; GPF77; GPR77; C5AR2
Species	Human
Gene ID (Human)	27202
Gene ID	27202
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 258-337 of human C5AR2 (NP_060955.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9P296
Protein Weight	36kDa
Shipping	Ice bag