

Anti-C5/C5A Polyclonal Antibody For Wb, Elisa - P01031

Item Number: CM0026256



Introduction

Rabbit polyclonal antibody targeting human C5/C5a, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for complement system and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	C5/C5a Rabbit pAb
Remarks/Alias	C5D; C5a; C5b; ECLZB; CPAMD4; C5/C5a
Species	Human
GeneID (Human)	727
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1567-1676 of human C5/C5a (NP_001726.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01031
Protein Weight	188kDa
Shipping	Ice bag