

Anti-C6Orf25 (G6B) Rabbit Polyclonal Antibody For Wb And Elisa - O95866

Item Number: CM0017542



Introduction

Rabbit polyclonal antibody against human C6orf25 (G6b), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for platelet and megakaryocyte inhibitory receptor studies.

[Learn More](#)

Parameter	Value
Product Name	C6orf25 Rabbit pAb
Remarks/Alias	G6b; NG31; G6b-B; THAMY; C6orf25
Species	Human
GeneID (Human)	80739
GeneID	80739
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-142 of human C6orf25 (NP_612121.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95866
Protein Weight	26kDa
Shipping	Ice bag