

Anti-C7 Polyclonal Antibody For Wb, If/Icc, Elisa - P10643

Item Number: CM0030387



Introduction

Rabbit polyclonal antibody targeting human C7, suitable for WB, IF/ICC, and ELISA. Detects human, mouse, and rat C7. Ideal for studying complement membrane attack complex formation and host defense mechanisms.

[Learn More](#)

Parameter	Value
Product Name	C7 Rabbit pAb
Remarks/Alias	C7
Species	Human
GeneID (Human)	730
GeneID	730
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 22-270 of human C7 (NP_000578.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10643
Protein Weight	94kDa
Shipping	Ice bag