

Anti-C7Orf50/Cholesin Polyclonal Antibody For Wb, Elisa - Q9Brj6

Item Number: CM0013701

Introduction

Rabbit polyclonal antibody against human C7orf50 (cholesin), raised against recombinant protein (aa 1-194). Validated for WB and ELISA; cross-reacts with mouse. Detects a band at ~22 kDa. Ideal for cholesterol metabolism and hormone signaling research.

[Learn More](#)

Parameter	Value
Product Name	C7orf50 Rabbit pAb
Remarks/Alias	YCR016W; C7orf50
Species	Human
GeneID	84310
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-194 of human C7orf50 (NP_115726.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BRJ6
Protein Weight	22kDa
Shipping	Ice bag