

Anti-C8Orf4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nr00

Item Number: CM0024891

Introduction

Rabbit polyclonal antibody targeting human C8orf4 (TC1/TC-1), suitable for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Regulates cell proliferation, differentiation, and inflammatory signaling.

[Learn More](#)

Parameter	Value
Product Name	C8orf4 Rabbit pAb
Remarks/Alias	TC1; TC-1; C8orf4
Species	Human
Gene ID (Human)	56892
Gene ID	56892
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-106 of human C8orf4 (NP_064515.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NR00
Protein Weight	12kDa
Shipping	Ice bag