


## Anti-C9 Rabbit Polyclonal Antibody For Wb, Elisa - P02748

Item Number: CM0030730



### Introduction

C9 Rabbit Polyclonal Antibody validated for Western blot and ELISA. Detects human and mouse complement component C9, a pore-forming MAC subunit. Suitable for innate immunity and complement pathway research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | C9 Rabbit pAb                                                                                                   |
| Remarks/Alias    | C9D; ARMD15; C9                                                                                                 |
| Species          | Human                                                                                                           |
| Gene ID (Human)  | 735                                                                                                             |
| Gene ID          | 735                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 22-265 of human C9 (NP_001728.1). |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, ELISA                                                                                                       |
| Cross Reactivity | Human, Mouse                                                                                                    |
| SWISS            | P02748                                                                                                          |
| Protein Weight   | 63kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |