


Anti-C9Orf114 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q5T280

Item Number: CM0021384



Introduction

Rabbit polyclonal antibody targeting human C9orf114 (CENP32/SPOUT1), validated for WB, IHC-P, ELISA. Cross-reacts with rat. Ideal for investigating centromere biology, mitotic spindle organization, and rRNA methylation.

[Learn More](#)

Parameter	Specification
Product Name	C9orf114 Rabbit pAb
Remarks/Alias	CENP32; CENP-32; HSPC109; C9orf114
Species	Human
Gene ID	51490 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 163-340 of human C9orf114 (NP_057474.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q5T280
Protein Weight	42kDa
Shipping	Ice bag