


Anti-C9orf72 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Lt7

Item Number: CM0023396



Introduction

Rabbit polyclonal antibody against human C9orf72 (aa 9-200). Suitable for WB, IHC-P, IF/ICC, and ELISA. Validated in human, mouse, and rat. Target involved in autophagy and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	C9orf72 Rabbit pAb
Remarks/Alias	ALSFTD, DENND9, FTDALS, DENNL72, FTDALS1, C9orf72
Species	Human
GeneID (Human)	203228
GeneID	203228
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 9-200 of human C9orf72 (NP_001242983.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96LT7
Protein Weight	54kDa
Shipping	Ice bag