


Anti-Ca13 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N1Q1

Item Number: CM0016589



Introduction

Rabbit polyclonal antibody against human carbonic anhydrase 13 (CA13), validated for WB, IHC-P, ELISA. Reacts with human and mouse. Essential tool for studying pH regulation and metabolic carbon dioxide conversion.

[Learn More](#)

Parameter	Value
Product Name	CA13 Rabbit pAb
Remarks/Alias	CAXIII; CA13
Species	Human
GeneID (Human)	377677
GeneID	377677
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human CA13 (NP_940986.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N1Q1
Protein Weight	29kDa
Shipping	Ice bag