


Anti-Ca3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P07451

Item Number: CM0017638



Introduction

High-quality rabbit polyclonal antibody targeting human carbonic anhydrase III (CA3). Validated for WB, IHC-P, and ELISA, with cross-reactivity to mouse and rat. Ideal for muscle-specific carbonic anhydrase detection. Weight ~30kDa.

[Learn More](#)

Parameter	Value
Product Name	CA3 Rabbit pAb
Remarks/Alias	Car3; CAIII; CA3
Species	Human
GeneID (Human)	761
GeneID	761
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human CA3 (NP_005172.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P07451
Protein Weight	30kDa
Shipping	Ice bag