

## Anti-Ca5A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P35218

Item Number: CM0017312

### Introduction

Rabbit polyclonal antibody targeting human CA5A (also mouse/rat). Validated for Western Blot, IF/ICC, and ELISA. Detects mitochondrial carbonic anhydrase crucial for urea cycle and Krebs cycle metabolism. For IVD research and diagnostic development.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | CA5A Rabbit pAb                                                                                                  |
| Remarks/Alias    | CA5; CAV; CAVA; CA5AD; GS1-21A4.1; CA5A                                                                          |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 763                                                                                                              |
| GeneID           | 763                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 39-305 of human CA5A (NP_001730.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IF/ICC, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | P35218                                                                                                           |
| Protein Weight   | 35kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |